

CLIK-FIX™

Luer Lock Securement

Clik-FIX Luer Lock Securement is designed to stabilize and secure any medical tubing junction with a luer lock connection. The foam adhesive base provides a cushion for the luer lock to minimize pressure points on the patient's skin. The foam top with luer lock window holds the luer lock in place by providing a physical barrier around the connection thereby reducing direct adhesive pressure on to the luer lock and patient.



- Center split release paper allows for easy removal of the release liner and placement of the device
- Skin friendly adhesive is solvent free, comfortable and holds well
- Stabilizes any medical tubing with luer lock connection
- Luer lock window feature reduces direct pressure on the device and patient
- Packaged sterile

ORDERING INFORMATION



LLS-1021

Clik-FIX Luer Lock Securement Device

50 per box



STARBOARD
MEDICAL

22805 - B Savi Ranch Parkway, Yorba Linda, CA 92887

P: 714.283.3099 F: 714.283.3033 www.starboardmedical.com

Starboard Medical™ and Clik-FIX™ are trademarks of Starboard Medical Inc, Yorba Linda, CA.

