


Comfortable, Soft, and Reliable

CLIK-FIX™ Foley Catheter Stabilization Device

Designed with "care in mind", Starboard Medical's Clik-FIX Foley Catheter Stabilization Device provides comfortable, soft, and reliable securement. The soft foam underpad conforms to the patient providing a cushion to reduce pressure points on the skin surface. The Wye Securement strap provides a solid barrier at the point where tugging most often occurs to help reduce catheter pistoning and accidental dislodgement.



Soft foam
underpad

Wye Securement
strap provides a solid
securement barrier

Velcro is adjustable
and repositionable

No hard
plastic parts

- Low profile securement
- No hard plastic parts
- Helps prevent urethral erosion, pistoning, and accidental dislodgement
- Adjustable, can be repositioned
- Easy to place and remove

ORDERING INFORMATION

Clik-FIX Foley Catheter Stabilization,
packaged 50 each per box.

FCS-1073	Secures 16/24 French
FCS-1074	Secures 10-14 French
FCS-1075	Secures 5-8 French



18001 Sky Park Circle, Suite E, Irvine, CA 92614

P: 714.283.3099 F: 714.283.3033 www.starboardmedical.com

Starboard Medical™ and Clik-FIX™ are trademarks of Starboard Medical Inc, Irvine, CA.

US Patent #9526869, other patents pending