


CLIK-FIX™

Neonatal PICC Securement Device



The Clik-FIX™ Neonatal PICC Securement device is designed with patient care in mind. It features a soft foam securement base with no hard plastic parts to cause pressure points, bruising or trauma to the baby's skin. The built in "Chevron" securement straps lay over the catheter wing to stabilize the catheter and minimize micropistoning and inadvertent dislodgement.

Clik-FIX Neonatal PICC Securement device features a universal design and is compatible with most neonatal PICC lines.

ORDERING INFORMATION

NPS-1006

Neonatal PICC Securement


50 per box

Features skin friendly hypoallergenic adhesive and center split release liner for easy and comfortable application


Soft comfortable foam cushions catheter wing reducing pressure ulcer risk

Velcro straps "Chevron" around the catheter wing to minimize catheter movement risk

Adhesive base gently holds catheter wing



It's as easy to place as one, two, three.



STARBOARD
MEDICAL

18001 Sky Park Circle, Suite E, Irvine, CA 92614

P: 714.283.3099 F: 714.283.3033 www.starboardmedical.com

Starboard Medical™ and Clik-FIX™ are trademarks of Starboard Medical Inc, Irvine, CA.

US Patent #9526869, other patents pending

