

CLIK-FIX™

Soft PICC / Central Securement Device



The Clik-FIX™ Soft PICC / Central Securement device is designed with patient care in mind. It features a soft foam securement base with no hard plastic parts to cause pressure points, bruising or trauma to the skin. The built in "Chevron" securement straps lay over the catheter wing to stabilize the catheter and minimize micropistoning and inadvertent dislodgement.

Clik-FIX Soft PICC / Central Securement device features a universal design and is compatible with most PICC and central lines.

ORDERING INFORMATION

WCS-1006

Soft PICC / Central Securement 50 per box


Features skin friendly hypoallergenic adhesive and center split release liner for easy and comfortable application

Soft comfortable foam cushions catheter wing reducing pressure ulcer risk

Velcro straps "Chevron" around the catheter wing to minimize catheter movement risk

Adhesive base gently holds catheter wing

It's as easy to place as one, two, three.



18001 Sky Park Circle, Suite E, Irvine, CA 92614

P: 714.283.3099 F: 714.283.3033 www.starboardmedical.com

Starboard Medical™ and Clik-FIX™ are trademarks of Starboard Medical Inc, Irvine, CA.

US Patent #9526869, other patents pending

