


# CLIK-FIX™

## All 'n One Peripheral Catheter Securement

Clik-FIX Peripheral Catheter Securement is an all-in-one Engineered Stabilization Device (ESD) with bio-occlusive dressing. Clik-FIX All 'n One Peripheral Catheter Securement replaces the need to use individual tape, gauze and transparent dressings. The window provides for a clear view of the insertion site and surrounding skin. Clik-FIX PCS is designed to cushion and stabilize the catheter hub and luer lock junction to help reduce the risk of device related pressure ulcers, catheter movement, dislodgement and associated complications. Superior securement methods improve outcomes and reduce costs.



- All 'n One Engineered Stabilization Device for Peripheral IV Catheters
- Bio-Occlusive dressing protects insertion site
- U-shaped window provides for visualization of access site and surrounding skin
- Foam base pad cushions the load on the skin at the catheter hub-luer lock junction reducing the risk of pressure injury
- Non-irritating adhesive foam and dressing provide superior securement yet are easy to remove
- Replaces tape and transparent dressing

### ORDERING INFORMATION



**PCS-1011**  
Clik-FIX All 'n One Peripheral IV Securement  
50 per box



18001 Sky Park Circle, Suite E, Irvine, CA 92614

P: 714.283.3099 F: 714.283.3033 [www.starboardmedical.com](http://www.starboardmedical.com)

Starboard Medical™ and Clik-FIX™ are trademarks of Starboard Medical Inc, Irvine, CA.

US Patent #9526869, other patents pending



