CLIK-FIX*

All 'n One Peripheral Catheter Securement

Clik-FIX Peripheral Catheter Securement is an all-in-one Engineered Stabilization Device (ESD) with bio-occlusive dressing. Clik-FIX All 'n One Peripheral Catheter Securement replaces the need to use individual tape, gauze and transparent dressings. The window provides for a clear view of the insertion site and surrounding skin. Clik-FIX PCS is designed to cushion and stabilize the catheter hub and luer lock junction to help reduce the risk of device related pressure ulcers, catheter movement, dislodgement and associated complications. Superior securement methods improve outcomes and reduce costs.



- All 'n One Engineered Stabilization Device for Peripheral IV Catheters
- · Bio-Occlusive dressing protects insertion site
- U-shaped window provides for visualization of access site and surrounding skin
- Foam base pad cushions the load on the skin at the catheter hub-luer lock junction reducing the risk of pressure injury
- Non-irritating adhesive foam and dressing provide superior securement yet are easy to remove
- · Replaces tape and transparent dressing

ORDERING INFORMATION



PCS-1011 Clik-FIX All 'n One Peripheral IV Securement 50 per box



18001 Sky Park Circle, Suite E, Irvine, CA 92614 P: 714.283.3099 F: 714.283.3033 www.starboardmedical.com



CLIK-FIX[™]

All 'n One Peripheral Catheter Securement

Designed to reduce the risk of pressure injury due to IV hubs

Features	Benefits	Clinical Advantages
All-in-one	Less bulk, less fussOne step dressing and catheter securement	 Replaces the need to use individual tape, gauze, transparent dressing, and a secondary securement device
U-shaped Window	 Provides a clear view of the insertion site and surrounding skin 	Bio-Occlusive dressing protects the insertion site from external contaminants
Luer Lock Holder	 Provides a barrier around bottom and top of luer lock to help prevent catheter movement 	 Stabilizes the catheter at the hub and also at the luer lock junction Provides two times the stabilization to help prevent catheter movement, dislodgement, and associated complications
Center Split Release Liner	 Allows for easy removal of the release liner and accurate placement of the device 	Easy to placeLess likely to stick to itself or gloves during application
Solvent Free Adhesive	Adhesive is not rough on skin	• Skin friendly and easy to remove • Improves patient comfort
Luer Lock Window	 Holds catheter hub and luer lock without placing direct pressure on junction 	 Helps prevent catheter movement and pistoning Reduces pressure applied on luer lock and hub
Reliable Adhesive Base	 Foam base pad cushions the load on the skin at the catheter hub-luer lock junction Holds catheter at slight angle 	 Non-irritating adhesive foam and dressing provide superior securement yet are easy to remove Helps maintain catheter insertion angle, prevent catheter movement, and reduce pressure ulcer risk



