


Tympanic Temperature Sensor

BARRIER FOAM, COLLAPSES AND EXPANDS SEALING OFF EAR CANAL FROM AMBIENT AIR PROVIDING ACCURATE READINGS WITHOUT PHYSICALLY TOUCHING TYMPANIC MEMBRANE



Tympanic Probe		
Catalog Number	Type	Quantity/Case
400-TY	Adult	50
400-TYP	Pediatric	50



18001 Sky Park Circle Suite E-F
Irvine, CA 92614

Phone: 714.283.3099
Fax: 714.283.3033

www.starboardmedical.com